

ABMAX premium



Premium Design - Premium Performance - Premium Results



Quick and precise volume setting


Adjusting the volume of the **LABMAX premium** is even simpler: just push the re-designed quick-adjust knob – slide to the desired volume– release – that's all!

Perfect air-purging mechanism for zero reagent loss

Turn the cannula at 90°. Remove (purge) trapped air without reagent loss by pumping in a closed circuit: the reagent flows back into the bottle!

After pumping 2 or 3 times the dispenser is free of air and ready to use.

This special feature not only helps to save residual reagent but also when accidentially setting the wrong volume!





Easy and fatigue-proof dispensing

Turn the cannula into the dosing mode and enjoy dispensing with **LABMAX premium**. The piston of the **LABMAX premium** completely discharges the reagent with each stroke. The dispensing glass cylinder is protected with a transparent plastic coating which allows the user to view the liquid in the glass cylinder. The protective coating also reduces

100% drip-free system for your safety

By turning the cannula at 180° all liquid in the discharge tube flows back into the bottle. The **LABMAX premium** is completely closed and no accidental discharge can occur by inadvertently moving the piston.

There is absolutely no dripping possible!





360°-turn for your added safety

the risk of breakage!

No matter how much you tighten the **LABMAX premium**, you can always redirect the cannula by simply turning the valve block to the desired position. This way the label of the bottle is always visible.

This feature works for all our adapters. Adapters for different bottle threads can stay on the bottle!

"Do-it-yourself" recalibration for a long life time

If you use the **LABMAX premium** as a measuring device in conformity to ISO 9001 and GLP, you may also want to control and monitor its accuracy and reliability. The piston can easily be recalibrated by turning an

internal screw!





To safe reagents, large volume glass instruments are replaced by small volume measuring instruments step by step. This reduces the costs and is a better alternative in terms of environmental aspects. **LABMAX premium** is our sophisticated dispensing system which fulfills your increasing requirements in laboratories:



Bottle-Top Dispenser Liquid Handling



ABMAX premium

Your advantages at a glance:

- all-purpose dispenser, universal chemical compatibility
- highest accuracy and reproducibility
- built-in recalibration mechanism
- robust and stable, made for a long life time
- zero loss of reagent by pumping in a closed circuit
- 100% drip-free due to our unique anti-drip system
- 360 degree rotatable valve block for added safety
- flexible by using different bottle adapters
- easy to clean, completely autoclavable at 121°C
- conformity certification and CE label
- Branded according to the German calibration law from 25.07.2013 page 2722-2723 on basis of the EU policy 2004/22/EG from 31.03.2004

Order information

LABMAX premium

Volume	Division	Tolerance	Deviation	Thread	Order number
ml	ml	ml	ml	GL	
0.25 - 2.5	0.05	±0.012	0.002	32	5 370 901 ⁼
0.5 - 5.0	0.10	±0.03	0.005	32	5 370 902
1.0 - 10.0	0.20	±0.06	0.010	32	5 370 903
2.5 - 25.0	0.50	±0.15	0.025	45	5 370 904
5.0 - 50.0	1.00	±0.30	0.050	45	5 370 905
10.0 - 100.0	2.00	±0.60	0.100	45	5 370 906
25.0 - 250.0	5.00	±1.60	0.250	45	5 370 907

Scope of supply

- filling tube
- discharge tube
- instruction manual
- conformity declaration
- quality certificate
- serial number

thread adapter 2.5/5.0/10.0ml: A25, A28, A32/45, S40

25.0/50.0/100.0ml: A25, A28, A45/32, S40

